

# PHYSICS PRACTICAL SHEETS

Date 09-03-2023

KV CAMPUS

Class CE

Roll No. 25

Shift Morning

Experiment No. 2

Group T

Sub. Physics

Set

Object of the Experiment (Block Letter)

DETERMINATION OF THE WAVELENGTH OF SODIUM LIGHT USING NEWTON'S RING APPARATUS

Apparatus Required:

- i) A travelling microscope
- ii) Sodium lamp
- iii) Newton's rings apparatus
- iv) A spherometer
- v) A convex lens of small focal length.

Theory:

Circular interference fringes produced by enclosing a thin air film of varying thickness between the surface of convex lens of large radius of curvature and a plane glass plate are known as Newton's rings. The wavelength of monochromatic light which produces these rings by

$$\lambda = \frac{D_n^2 - D_m^2}{4R(n-m)}$$

where,  $R$  = radius of curvature of surface of plane convex lens in contact with the glass plate.

$D_n$  and  $D_m$  are diameters of  $n^{\text{th}}$  and  $m^{\text{th}}$  dark/bright rings.

Observations:

Least count of main scale =  $0.1 \text{ cm}$

No. of vernier scale division =  $100$

Vernier constant (V.C) =  $\frac{0.1}{100} = 0.001 \text{ cm}$




Fig: Schematic diagram of  
Newton's ring apparatus



← D →

Fig: Newton's rings.

Table for diameter of the Newton's ring

No. of obs	Ring No	Microscope reading		Diameter	Microscope Read.		Diameter	Mean
		Left	Right		Left	Right		
1	20	4.330	3.330	1	3.335	4.328	0.993	0.99
2	16	4.279	3.040	0.874	3.385	4.280	0.895	0.887
3	12	4.222	3.45	0.712	3.442	4.222	0.780	0.776
4	8	4.154	3.52	0.634	3.511	4.155	0.644	0.639
5	4	4.069	3.61	0.459	3.591	4.070	0.479	0.469

Result:

The radius of curvature  $R$  of plano-convex lens ( $R$ ) =  $\frac{a^2}{6h} + \frac{h}{2}$   
 $= 260 \text{ cm}$

The wavelength  $\lambda$  obtained from  $D_{20}, D_2 = 4.709 \times 10^{-5} \text{ cm}$

The wavelength  $\lambda$  obtained from  $D_{16}, D_8 = 4.643 \times 10^{-5} \text{ cm}$

The wavelength  $\lambda$  obtained from  $D_{12}, D_4 = 4.649 \times 10^{-5} \text{ cm}$

The wavelength of sodium light is found to be  
 $= 4.649 \times 10^{-5} \text{ cm}$

$$\% \text{ error} = \left| \frac{5.8 \times 10^{-5} - 4.649 \times 10^{-5}}{5.8 \times 10^{-5}} \right| \times 100\%$$

$$= 19.84\%$$

### Precautions:

- (i) The glass plate and lens must be cleaned properly.
- (ii) The lens must have larger radius of curvature.
- (iii) The amount of light from the source should be adjusted for maximum visibility.